

7 M under

At main Buss which is 90' away from location of ~~the~~ instrument

- ① 208/110 volts with no measurable fluctuations
- ② No data on 30 Amp load drop.
- ③ Measured no load voltage at instrument site 206+/109-volts

I. Permutation Phoned to

STAT

STAT

STAT

INSTALLATION ENGINEERINGFile
PI-624

I. INSTRUMENT

A. Name Enlarger, Coherent

B. Manufacturer _____

C. Contract Number _____

STAT

II. PHYSICAL FEATURES

A. Number of Component Parts Three (enlarger, elec. rack, compressor)

B. Dimensions of the Largest Component Part:

Length 11 Ft. - In. Height 6 Ft. - In.Width 3 Ft. 10 In. NOTE: Height of equip. rack = 6'10"C. Weight of Largest Component Part 2,150 lbs.D. Total Weight of Instrument 2,750

E. Overall Dimensions Assembled:

Length _____ Ft. _____ In. Does not Apply Height _____ Ft. _____ In.

Width _____ Ft. _____ In.

F. Type of Base of Mount:

Flat X Three Point Suspension _____ Four Point Suspension _____G. Does Instrument have built-in mobility? NoH. Is the instrument particularly sensitive to vibration? No

I. Are any special or unusual tools or fixtures necessary or advisable for the installation or maintenance of this equipment? _____

Installation will be accomplished by Perkin-Elmer personnel & equipment.

III. UTILITIES

A. Electrical:

Voltage

AC
110 Volts + 5 Volts

DC

None

Current

30 Amps Max.

Frequency

60 cps

Nr. of phases

one

Nr. of wires

three (incl. ground)

Power required by equipment

3K Wattsnone Watts

Type of outlet required: Two Prong _____, Three Prong _____

Twist Lock _____, Permanent Installation X

Should the equipment be shielded, either from external electro-magnetic signals, or to prevent interference with other equipment?

No

B. Air Conditioning:

Room temperature 70°F + 10°F Humidity No control required
 Output of Instrument 10,000 BTU/Hr. incl. air compressor
 If air must be filtered, what is maximum permissible particle size
 in microns? See reverse side What particle count? _____
 particles per cubic foot.
 Direct connection to instrument? Yes _____ No X
 If yes to above, what is the desired air temperature to instrument? _____
 Should discharged air be ducted separately? Yes
 Is discharged air noxious? Yes toxic? Yes
 Connector size to instrument 4" dia.

C. Plumbing:

Is water required for the instrument? Yes _____ No X
 Water pressure _____ Flow in GPM _____
 Type of water desired:
 Tap _____ °F + _____ °F
 Tempered _____ °F + _____ °F
 Deionized _____ °F + _____ °F
 Filtered _____ °F + _____ °F Particle size and count per
 unit volume.
 Type of pipe required:
 Galvanized _____ Copper _____
 Stainless Steel _____ Plastic _____
 Is floor drain required? Yes _____ No _____
 Diameter of drain _____ Galvanized drain _____
 Plastic drain _____ Glass drain _____

D. Compressed Air: See reverse side.

Diameter of connectors _____ Type of connectors _____
 PSI _____ Water free? _____
 CFM _____ Oil free? _____

E. Vacuum: See reverse side.

Is vacuum required? Yes _____ No _____
 Vacuum required _____ PSIA or _____ (inches) (milli-
 meters) of Hg
 Displacement _____ CFM _____

IV. REMARKS

In the event additional space is required for environmental conditions
 or utilities not mentioned above, use the reverse side of this form.

Air conditioning: Clean room operation is recommended as dust on lenses is detrimental. Operator can determine exact requirements by experience.

Note: Air Compressor and Vacuum Pump will be supplied with the instrument. Air Compressor will include sufficient hose and cable to permit location in adjacent room. Air compressor has automatic water drain, will drain into catch bucket.